


Mbd2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ubb5

Item Number: CM0017224



Introduction

Rabbit polyclonal antibody against human MBD2, a methyl-CpG-binding transcription repressor and NuRD complex component. Validated for WB and ELISA, cross-reacts with human, mouse, rat. Ideal for epigenetic and gene silencing studies.

[Learn More](#)

Parameter	Value
Product Name	MBD2 Rabbit pAb
Remarks/Alias	DMTase; NY-CO-41; MBD2
Species	Human
GeneID	8932
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-411 of human MBD2 (NP_003918.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBB5
Protein Weight	43kDa
Shipping	Ice bag