

Mbd6 Rabbit Pab For Wb, Elisa - Q96Dn6

Item Number: CM0024213

Introduction

Anti-MBD6 rabbit polyclonal antibody validated for WB and ELISA, reacting with human, mouse, and rat. Suitable for epigenetics and chromatin research applications.

[Learn More](#)

Parameter	Value
Product Name	MBD6 Rabbit pAb
Remarks/Alias	MBD6
Species	Human
GeneID (Human)	114785
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 880-980 of human MBD6 (NP_443129.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96DN6
Protein Weight	101kDa
Shipping	Ice bag