

Mcf2 Rabbit Polyclonal Antibody For Wb And Elisa - P10911

Item Number: CM0014656

Introduction

Rabbit polyclonal antibody targeting human MCF2 (proto-oncogene DBL), suitable for Western blot and ELISA applications. Cross-reacts with mouse. Useful for studying Rho GTPase signaling and oncogenic pathways.

[Learn More](#)

Parameter	Value
Product Name	MCF2 Rabbit pAb
Remarks/Alias	DBL; ARHGEF21; MCF2
Species	Human
Gene ID (Human)	4168
Gene ID	4168
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 872-1001 of human MCF2 (NP_001165347.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10911
Protein Weight	108kDa
Shipping	Ice bag