

Mcl1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q07820

Item Number: CM0016067



Introduction

High-quality rabbit polyclonal antibody targeting human MCL1 for Western blot, immunofluorescence, and ELISA applications. Cross-reacts with mouse and rat. Suitable for apoptosis and cell survival studies.

[Learn More](#)

Parameter	Value
Product Name	MCL1 Rabbit pAb
Remarks/Alias	TM; EAT; MCL1L; MCL1S; Mcl-1; BCL2L3; MCL1-ES; bcl2-L-3; mcl1/EAT; MCL1
Species	Human
GeneID (Human)	4170
GeneID	4170
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q07820
Protein Weight	37kDa
Shipping	Ice bag