

Mcl1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q07820

Item Number: CM0013406



Introduction

MCL1 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human and mouse MCL1 (37 kDa), a key apoptosis regulator. Ideal for cell survival and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MCL1 Rabbit pAb
Remarks/Alias	TM; EAT; MCL1L; MCL1S; Mcl-1; BCL2L3; MCL1-ES; bcl2-L-3; mcl1/EAT; MCL1
Species	Human
GeneID (Human)	4170
GeneID	4170
Immunogen	Synthetic peptide within amino acids 150-250 of human MCL1 (NP_068779.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q07820
Protein Weight	37kDa
Shipping	Ice bag