


Mcm4 Rabbit Polyclonal Antibody For Wb And Elisa - P33991

Item Number: CM0018054



Introduction

Rabbit polyclonal antibody targeting human MCM4, validated for WB and ELISA with reactivity to human, mouse, and rat. Essential tool for studying DNA replication licensing and cell cycle regulation.

[Learn More](#)

Parameter	Value
Product Name	MCM4 Rabbit pAb
Remarks/Alias	NKCD; CDC21; CDC54; IMD54; NKGCD; hCdc21; P1-CDC21; MCM4
Species	Human
Gene ID (Human)	4173
Gene ID	4173
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33991
Protein Weight	97kDa
Shipping	Ice bag