


Mcm5 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - P33992

Item Number: CM0012622



Introduction

Rabbit monoclonal antibody targeting human MCM5, a DNA replication licensing factor. Validated for WB, IHC-P, IP, and ELISA applications, with cross-reactivity to mouse and rat. Suitable for cell cycle and DNA replication research.

[Learn More](#)

Parameter	Value
Product Name	MCM5 Rabbit mAb
Remarks/Alias	CDC46; MGORS8; P1-CDC46; MCM5
Species	Human
GeneID (human)	4174
GeneID	4174
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 635-734 of human MCM5 (P33992).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33992
Protein Weight	82kDa
Shipping	Ice bag