


Mcm7 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P33993

Item Number: CM0015890



Introduction

MCM7 rabbit polyclonal antibody, validated for western blot, immunofluorescence, immunoprecipitation, and ELISA. Reacts with human, mouse, rat. Suitable for DNA replication and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	MCM7 Rabbit pAb
Remarks/Alias	MCM2; CDC47; P85MCM; P1CDC47; PNAS146; PPP1R104; P1.1-MCM3
Species	Human
Gene ID (Human)	4176
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human MCM7 (NP_005907.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33993
Protein Weight	81kDa
Shipping	Ice bag