


Mda5 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Byx4

Item Number: CM0004605



Introduction

Rabbit monoclonal antibody targeting human MDA5 (IFIH1), a cytoplasmic sensor of viral RNA. Validated for WB and ELISA, it cross-reacts with human, mouse, and rat. Ideal for innate immunity and infectious disease research.

[Learn More](#)

Parameter	Specification
Product Name	MDA5 Rabbit mAb
Remarks/Alias	AGS7; Hlcd; MDA5; IMD95; MDA-5; RLR-2; IDDM19; SGMRT1
Species	Human
GeneID (Human)	64135
GeneID	64135
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 900-1000 of human MDA5 (Q9BYX4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYX4
Protein Weight	25 kDa/117 kDa
Shipping	Ice bag