


Mda5 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q9Byx4

Item Number: CM0009083



Introduction

MDA5 Rabbit Monoclonal Antibody (Q9BYX4) for WB, IP, ELISA. Detects human IFIH1, a cytoplasmic sensor of viral nucleic acids. Suitable for innate immune response and viral infection research.

[Learn More](#)

Parameter	Value
Product Name	MDA5 Rabbit mAb
Remarks/Alias	AGS7; Hlcd; MDA5; IMD95; MDA-5; RLR-2; IDDM19; SGMRT1
Species	Human
GeneID	64135
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human IFIH1 (NP_071451.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q9BYX4
Protein Weight	25 kDa / 117 kDa
Shipping	Ice bag