


Mda5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Byx4

Item Number: CM0017746



Introduction

Rabbit polyclonal antibody against MDA5 (IFIH1), a cytoplasmic viral RNA sensor. Validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Ideal for innate immunity and antiviral research.

[Learn More](#)

Parameter	Value
Product Name	MDA5 Rabbit pAb
Remarks/Alias	AGS7; Hlcd; MDA5; IMD95; MDA-5; RLR-2; IDDM19; SGMRT1; IFIH1
Species	Human
GeneID (Human)	64135
GeneID	64135
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human IFIH1 (NP_071451.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYX4
Protein Weight	25 kDa/117 kDa
Shipping	Ice bag