

Mdfi Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q99750

Item Number: CM0017577

Introduction

MDFI Rabbit pAb validated for Western blot, IF/ICC, and ELISA applications. Detects endogenous MDFI in human, mouse, and rat samples. MyoD family inhibitor involved in myogenesis and Wnt signaling.

[Learn More](#)

Parameter	Value
Product Name	MDFI Rabbit pAb
Remarks/Alias	I-MF; I-mfa; MDFI
Species	Human
GeneID (Human)	4188
GeneID	4188
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human MDFI (NP_005577.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99750
Protein Weight	25kDa
Shipping	Ice bag