

Mdh2 Rabbit Pab - P40926

Item Number: CM0014570



Introduction

High-quality anti-MDH2 polyclonal antibody for WB and ELISA. Specific for human and mouse mitochondrial malate dehydrogenase, essential for TCA cycle and malate-aspartate shuttle studies. Ideal for metabolic and mitochondrial disorder research.

[Learn More](#)

Parameter	Value
Product Name	MDH2 Rabbit pAb
Remarks/Alias	MDH, MOR1, DEE51, M-MDH, EIEE51, MGC:3559, MDH2
Species	Human
Gene ID (Human)	4191
Gene ID	4191
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 59-338 of human MDH2 (NP_005909.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P40926
Protein Weight	36kDa
Shipping	Ice bag