

Med25 Rabbit Pab - Q71Sy5

Item Number: CM0019723



Introduction

Rabbit polyclonal antibody targeting human MED25 (ARC92/p78), validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for investigating Mediator complex and transcriptional regulation.

[Learn More](#)

Parameter	Specification
Product Name	MED25 Rabbit pAb
Remarks/Alias	P78; ACID1; ARC92; BVSYS; PTOV2; CMT2B2; TCBAP0758; MED25
Species	Human
Gene ID (Human)	81857
Gene ID	81857
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 30-200 of human MED25 (NP_112235.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q71SY5
Protein Weight	78kDa
Shipping	Ice bag