

Med26 Rabbit Pab For Wb And Elisa - O95402

Item Number: CM0020943



Introduction

MED26 Rabbit pAb is a polyclonal antibody raised in rabbit against recombinant human MED26 protein. It detects MED26 in human, mouse, and rat by WB and ELISA, aiding transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	MED26 Rabbit pAb
Remarks/Alias	CRSP7; CRSP70; MED26
Species	Human
GeneID (Human)	9441
GeneID	9441
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 500-574 of human MED26 (NP_004822.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95402
Protein Weight	65kDa
Shipping	Ice bag