

Mef2A+Mef2C Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q06413 / Q02078

Item Number: CM0012635

Introduction

Rabbit monoclonal antibody targeting MEF2A and MEF2C, validated for WB, IHC-P, and ELISA in human and mouse. Ideal for studying muscle differentiation, neuronal development, and immune cell function.

[Learn More](#)

Parameter	Value
Product Name	MEF2A+MEF2C Rabbit mAb
Remarks/Alias	NEDHSIL; DEL5q14.3; C5DELq14.3; mef2; ADCAD1; RSRFC4; RSRFC9; MEF2A+MEF2C
Species	Human
GeneID (Human)	4208/4205
GeneID	4208/4205
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 250-473 of human MEF2A+MEF2C (Q06413).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q06413/Q02078
Protein Weight	51kDa
Shipping	Ice bag