

## Megf10 Rabbit Pab - Q96Kg7

Item Number: CM0015482

### Introduction

MEGF10 Rabbit pAb is a polyclonal antibody targeting human MEGF10 (Q96KG7). Validated for Western blot and ELISA, it detects endogenous MEGF10 in human, mouse, and rat samples. Ideal for phagocytosis and muscle biology research.

[Learn More](#)

| Parameter        | Specification                                                                                                      |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | MEGF10 Rabbit pAb                                                                                                  |
| Remarks/Alias    | SR-F3, EMARDD, CMYP10A, CMYP10B, MEGF10                                                                            |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 84466                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 26-160 of human MEGF10 (NP_115822.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q96KG7                                                                                                             |
| Protein Weight   | 122kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                            |