

Mek3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P46734

Item Number: CM0025852

Introduction

High-quality MEK3 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, and ELISA applications. Reacts with human, mouse, and rat MEK3 (MAP2K3). Ideal for stress signaling and MAPK pathway research.

[Learn More](#)

Parameter	Value
Product Name	MEK3 Rabbit pAb
Remarks/Alias	MEK3; MKK3; MAPKK3; PRKMK3; SAPKK2; SAPKK-2
Species	Human
GeneID (Human)	5606
GeneID	5606
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human MEK3 (NP_659731.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46734
Protein Weight	39kDa
Shipping	Ice bag