

Mek3/Mek6 Rabbit Monoclonal Antibody For Wb And Elisa - P46734 / P52564

Item Number: CM0003792

Introduction

Rabbit monoclonal antibody targeting human MEK3/MEK6. Validated for WB and ELISA, cross-reacts with mouse and rat. Designed for studying MAP kinase signaling pathways in stress and cytokine responses.

[Learn More](#)

Parameter	Value
Product Name	MEK3/MEK6 Rabbit mAb
Remarks/Alias	MAPKK3; MEK3; MKK3; PRKMK3; SAPKK-2; SAPKK2; MEK3/MEK6
Species	Human
GeneID (Human)	5606/5608
GeneID	5606/5608
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 248-347 of human MEK3/MEK6 (P46734).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46734/P52564
Protein Weight	39,37kDa
Shipping	Ice bag