


Menin Rabbit Monoclonal Antibody For Wb, Elisa - O00255

Item Number: CM0010504



Introduction

High-quality rabbit monoclonal Menin antibody validated for Western blot and ELISA. Reacts with human, mouse, and rat Menin. Suitable for chromatin biology and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	Menin Rabbit mAb
Remarks/Alias	MEAI; SCG2; Menin
Species	Human
GeneID (human)	4221
GeneID	4221
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 516-615 of human Menin (O00255).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00255
Protein Weight	68kDa
Shipping	Ice bag