

Merlin Rabbit Polyclonal Antibody For Wb, Elisa - P35240

Item Number: CM0022620



Introduction

Rabbit polyclonal antibody against human Merlin (NF2), validated for WB and ELISA. Detects endogenous Merlin in human, mouse, and rat samples. Suitable for Hippo signaling, tumor suppression, and membrane biology studies.

[Learn More](#)

Parameter	Value
Product Name	Merlin Rabbit pAb
Remarks/Alias	ACN; SCH; BANF; merlin-1; Merlin
Species	Human
GeneID (Human)	4771
GeneID	4771
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 347-595 of human Merlin (NP_000259.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35240
Protein Weight	70kDa
Shipping	Ice bag