

Metap2 Rabbit Monoclonal Antibody For Wb, Elisa - P50579

Item Number: CM0008207



Introduction

METAP2 Rabbit Monoclonal Antibody validated for WB and ELISA. Recognizes human, mouse, and rat METAP2, a methionine aminopeptidase critical for N-terminal methionine removal and regulation of protein synthesis.

[Learn More](#)

Parameter	Value
Product Name	METAP2 Rabbit mAb
Remarks/Alias	MAP2; MNPEP; p67eIF2; METAP2
Species	Human
GeneID (Human)	10988
GeneID	10988
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 380-478 of human METAP2 (NP_006829.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50579
Protein Weight	53kDa
Shipping	Ice bag