

Methyl-Hp1 Gamma/Cbx3 (H175) Rabbit Polyclonal Antibody Q13185 For Dot Blot & Elisa

Item Number: CM0000222



Introduction

Rabbit polyclonal antibody targeting methylated HP1 gamma/CBX3 at H175, validated for dot blot and ELISA in human samples. Recognizes human CBX3 protein (21 kDa). Suitable for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Methyl-HP1 gamma/CBX3-H175 Rabbit pAb
Remarks/Alias	HECH; HP1-GAMMA; HP1Hs-gamma; Methyl-HP1 gamma/CBX3-H175
Species	Human
GeneID (human)	11335
GeneID	11335
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-183 of human HP1 gamma/CBX3 (NP_009207.2).
Source	Rabbit
Category	Methylated Antibodies
Application	DB, ELISA
Cross Reactivity	Human
SWISS	Q13185
Protein Weight	21kDa
Shipping	Ice bag