

Metrnl Rabbit Polyclonal Antibody For Wb, Elisa - Q641Q3

Item Number: CM0021953

Introduction

Rabbit polyclonal antibody targeting Meteorin-like protein (METRNL). Validated for Western blot (WB) and ELISA applications. Recognizes human and mouse METRNL. Ideal for research in energy expenditure, adipose tissue browning, and metabolic regulation.

[Learn More](#)

Parameter	Value
Product Name	METRNL Rabbit pAb
Species	Mouse
Gene ID	284207
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 179-311 of mouse METRNL. (NP_001004431.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q641Q3
Protein Weight	34kDa
Shipping	Ice bag