

Mettl14 Rabbit Pab - Q9Hce5

Item Number: CM0025238



Introduction

Rabbit anti-METTL14 polyclonal antibody, key component of m6A RNA methylation complex. Validated for WB, IF-P, IHC-P, ELISA in human, mouse, rat. Ideal for epitranscriptomics and RNA biology research.

[Learn More](#)

Parameter	Value
Product Name	METTL14 Rabbit pAb
Remarks/Alias	hMETTL14; METTL14
Species	Human
GeneID (Human)	57721
GeneID	57721
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human METTL14 (NP_066012.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HCE5
Protein Weight	52kDa
Shipping	Ice bag