


Mettl16 Rabbit Pab For Wb, Ip, If/Icc, Ihc-P, Elisa - Q86W50

Item Number: CM0016069



Introduction

High-quality rabbit polyclonal antibody targeting human METTL16 (RNA N6-methyltransferase). Validated for WB, IP, IF/ICC, IHC-P, and ELISA. Ideal for studying RNA modification and S-adenosyl-L-methionine homeostasis.

[Learn More](#)

Parameter	Value
Product Name	METTL16 Rabbit pAb
Remarks/Alias	METT10D; METTL16
Species	Human
GeneID (Human)	79066
GeneID	79066
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 263-562 of human METTL16 (NP_076991.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q86W50
Protein Weight	64 kDa
Shipping	Ice bag