

Mettl5 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nrn9

Item Number: CM0025502



Introduction

METTL5 Rabbit Polyclonal Antibody validated for WB, ELISA in human, mouse, rat. Targets rRNA N6-adenosine-methyltransferase METTL5 (Q9NRN9), essential for translation and stem cell pluripotency research.

[Learn More](#)

Parameter	Value
Product Name	METTL5 Rabbit pAb
Remarks/Alias	MRT72; HSPC133; METTL5
Species	Human
Gene ID (Human)	29081
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-209 of human METTL5 (NP_054887.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRN9
Protein Weight	24kDa
Shipping	Ice bag