

MgII Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - Q99685

Item Number: CM0007120



Introduction

High-quality rabbit monoclonal antibody against human MGLL (monoglyceride lipase, MAGL) for IVD research, validated for WB, IF-P, IHC-P, and ELISA, with cross-reactivity to mouse and rat. Ideal for endocannabinoid signaling and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	MGLL Rabbit PolymAb®
Remarks/Alias	MGL; HUK5; MAGL; HU-K5
Species	Human
Gene ID (Human)	11343
Gene ID	11343
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-303 of human MGLL. (NP_001003794.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99685
Protein Weight	33kDa
Shipping	Ice bag