


Mglur2 Rabbit Monoclonal Antibody For Wb - Q14416

Item Number: CM0009357



Introduction

mGluR2 Rabbit Monoclonal Antibody (Q14416) validated for WB, IF-P, and ELISA. Detects human metabotropic glutamate receptor 2, with cross-reactivity to mouse and rat. Ideal for neuroscience research on synaptic transmission and neuronal excitability.

[Learn More](#)

Parameter	Value
Product Name	mGluR2 Rabbit mAb
Remarks/Alias	GLUR2; mGluR2
Species	Human
Gene ID (Human)	2912
Gene ID	2912
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 773-872 of human mGluR2 (Q14416).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14416
Protein Weight	96 kDa
Shipping	Ice bag