


Mgmt Rabbit Polyclonal Antibody For Wb And Elisa - P16455

Item Number: CM0012817



Introduction

Rabbit polyclonal antibody targeting human MGMT (methylated-DNA--protein-cysteine methyltransferase). Validated for WB and ELISA applications. Useful for DNA repair and drug resistance studies.

[Learn More](#)

Parameter	Value
Product Name	MGMT Rabbit pAb
Remarks/Alias	MGMT
Species	Human
GeneID (human)	4255
GeneID	4255
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-207 of human MGMT (NP_002403.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P16455
Protein Weight	22kDa
Shipping	Ice bag