

Mhc Class I (Hla-A/B) Rabbit Mab -
P30443/P01889/P01891/P01892/P13746/P16188/P30450

Item Number: CM0008170

Introduction

Recombinant rabbit monoclonal antibody targeting human MHC class I (HLA-A/B). Validated for WB, IF/ICC, and ELISA. Suitable for immunology research and diagnostic assay development targeting HLA alleles.

[Learn More](#)

Parameter	Value
Product Name	MHC class I (HLA-A/B) Rabbit mAb
Remarks/Alias	HLAA; MHC class I (HLA-A/B)
Species	Human
Gene ID (Human)	3105/3106
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human MHC class I (HLA-A/B) (P04439).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P30443/P01889/P01891/P01892/P13746/P16188/P30450
Protein Weight	41kDa
Shipping	Ice bag