


Mib1/Dip1 Rabbit Monoclonal Antibody For Wb, Elisa - Q86Yt6

Item Number: CM0012549



Introduction

High-specificity rabbit monoclonal antibody targeting MIB1/DIP1 (Q86YT6). Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Suitable for studying Notch signaling, ubiquitination, and centrosome biology.

[Learn More](#)

Parameter	Value
Product Name	MIB1/DIP1 Rabbit mAb
Remarks/Alias	MIB; DIP1; ZZZ6; DIP-1; LVNC7; ZZANK2; MIB1/DIP1
Species	Human
GeneID (Human)	57534
GeneID	57534
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human MIB1/DIP1 (Q86YT6)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86YT6
Protein Weight	110kDa
Shipping	Ice bag