

Mica Rabbit Mab For Wb, Elisa - Q29983

Item Number: CM0003447



Introduction

Rabbit monoclonal antibody against human MICA (MHC class I polypeptide-related sequence A). Suitable for Western blot and ELISA; cross-reacts with human and mouse. Research use in immune surveillance and stress response studies.

[Learn More](#)

Parameter	Value
Product Name	MICA Rabbit mAb
Remarks/Alias	MIC-A; PERB11.1; MICA
Species	Human
GeneID (Human)	100507436
GeneID	100507436
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q29983
Protein Weight	43kDa
Shipping	Ice bag