

Mica Rabbit Pab For Wb, Elisa - Q29983

Item Number: CM0016251



Introduction

Polyclonal MICA antibody validated for WB and ELISA. Reacts with human, mouse, and rat. Targets MHC class I-related protein involved in NK cell activation; ideal for immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	MICA Rabbit pAb
Remarks/Alias	MIC-A; PERB11.1; MICA
Species	Human
GeneID	100507436
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-307 of human MICA (NP_000238.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q29983
Protein Weight	43kDa
Shipping	Ice bag