

Mical1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Tdz2

Item Number: CM0015829



Introduction

MICAL1 Rabbit pAb is a polyclonal antibody validated for WB and ELISA, targeting human and mouse MICAL1, a monooxygenase regulating actin dynamics and apoptosis. Suitable for research in cytoskeleton and neurobiology.

[Learn More](#)

Parameter	Value
Product Name	MICAL1 Rabbit pAb
Remarks/Alias	MICAL; NICAL; MICAL-1; MICAL1
Species	Human
Gene ID (Human)	64780
Gene ID	64780
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 718-1067 of human MICAL1 (NP_073602.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8TDZ2
Protein Weight	118kDa
Shipping	Ice bag