


Mitofusin 2 Rabbit Mab For Wb - O95140

Item Number: CM0002600



Introduction

Mitofusin 2 Rabbit mAb (O95140) validated for WB, IF-P, and ELISA, with cross-reactivity to human, mouse, and rat. A critical reagent for studying mitochondrial dynamics, fusion, and mitophagy.

[Learn More](#)

Parameter	Value
Product Name	Mitofusin 2 Rabbit mAb
Remarks/Alias	HSG; MARF; CMT2A; CPRP1; CMT2A2; HMSN6A; CMT2A2A; CMT2A2B; Mitofusin 2
Species	Human
GeneID (Human)	9927
GeneID	9927
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Mitofusin 2 (O95140).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95140
Protein Weight	86kDa
Shipping	Ice bag