

Mitofusin 2 Rabbit Pab - Mitofusin 2

Item Number: CM0012830

Introduction

Mitofusin 2 Rabbit Polyclonal Antibody for WB, IHC-P, ELISA applications. Detects Human, Mouse, Rat. Recombinant fusion protein immunogen corresponding to amino acids 70-260 of mouse Mitofusin 2.

[Learn More](#)

Parameter	Value
Product Name	Mitofusin 2 Rabbit pAb
Remarks/Alias	Fzo; D630023P19Rik; Mitofusin 2
Species	Mouse
Gene ID (Human)	9927
Gene ID	170731
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-260 of mouse Mitofusin 2 (NP_573464.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
Protein Weight	68kDa/86kDa
Shipping	Ice bag