

MLH3 Rabbit Pab For Wb And Elisa - Q9Uhc1

Item Number: CM0030007

Introduction

MLH3 Rabbit pAb is a polyclonal antibody targeting the DNA mismatch repair protein MLH3. Suitable for Western blot and ELISA in human, mouse, and rat samples, it supports research into DNA repair mechanisms and cancer.

[Learn More](#)

Parameter	Value
Product Name	MLH3 Rabbit pAb
Remarks/Alias	HNPCC7; MLH3
Species	Human
GeneID	27030
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1190-1389 of human MLH3 (NP_055196.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHC1
Protein Weight	164kDa
Shipping	Ice bag