


Mmp13 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P45452

Item Number: CM0013566



Introduction

MMP13 Rabbit Polyclonal Antibody (P45452) validated for Western blot, IF/ICC, and ELISA in human, mouse, and rat samples. Targets collagenase 3 essential for ECM degradation, bone development, and tissue remodeling.

[Learn More](#)

Parameter	Value
Product Name	MMP13 Rabbit pAb
Remarks/Alias	CLG3; MDST; MANDP1; MMP-13; 13
Species	Human
GeneID (human)	4322
GeneID	4322
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 262-471 of human MMP13 (NP_002418.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45452
Protein Weight	54kDa
Shipping	Ice bag