

Mmp17 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Ulz9

Item Number: CM0029033

Introduction

Highly specific MMP17 rabbit polyclonal antibody validated for WB, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for studying extracellular matrix remodeling and tumor biology.

[Learn More](#)

Parameter	Value
Product Name	MMP17 Rabbit pAb
Remarks/Alias	MMP-17; MT4MMP; MTMMP4; MT4-MMP; MMP17
Species	Human
GeneID (human)	4326
GeneID	4326
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9ULZ9
Protein Weight	67kDa
Shipping	Ice bag