

Mmp19 Rabbit Polyclonal Antibody For Wb, Elisa - Q99542

Item Number: CM0029106

Introduction

Rabbit polyclonal antibody targeting MMP19 (Matrix metalloproteinase-19). Validated for WB, ELISA. Detects human, mouse, rat MMP19; suitable for ECM degradation and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	MMP19 Rabbit pAb
Remarks/Alias	CODA; MMP18; RASI-1; MMP19
Species	Human
Gene ID (Human)	4327
Gene ID	4327
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 209-508 of human MMP19 (NP_002420.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q99542
Protein Weight	57kDa
Shipping	Ice bag