


Mmp9 Rabbit Polyclonal Antibody For Wb, Elisa - P14780

Item Number: CM0022770



Introduction

High-quality rabbit polyclonal antibody against human MMP9 (Gelatinase B). Validated for WB and ELISA, cross-reacts with mouse MMP9. Ideal for extracellular matrix research and IVD development.

[Learn More](#)

Parameter	Value
Product Name	MMP9 Rabbit pAb
Remarks/Alias	GELB; CLG4B; MMP-9; MANDP2; MMP9
Species	Human
Gene ID (Human)	4318
Gene ID	4318
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 538-707 of human MMP9 (NP_004985.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P14780
Protein Weight	78 kDa
Shipping	Ice bag