

Mneongreen Rabbit Polyclonal Antibody - Mneongreen

Item Number: CM0020249



Introduction

Rabbit polyclonal antibody against mNeonGreen fluorescent protein (Branchiostoma lanceolatum). Suitable for WB, IF/ICC, and ELISA applications with species-independent cross-reactivity. Expected band size ~30 kDa. For research use only.

[Learn More](#)

Parameter	Value
Product Name	mNeonGreen Rabbit pAb
Species	Branchiostoma lanceolatum
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-236 of Branchiostoma lanceolatum mNeonGreen.
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Species independent
Protein Weight	30kDa
Shipping	Ice bag