

Mnx1/Hb9/Hlxb9 Rabbit Monoclonal Antibody For Wb And Elisa - P50219

Item Number: CM0003717

Introduction

Rabbit monoclonal antibody against MNX1/HB9/HLXB9, raised with synthetic peptide. Validated for WB and ELISA with human and mouse samples. Suitable for studying motor neuron identity and pancreas development.

[Learn More](#)

Parameter	Value
Product Name	MNX1/HB9/HLXB9 Rabbit mAb
Remarks/Alias	HB9; HLXB9; SCRA1; HOXHB9; MNX1/HB9/HLXB9
Species	Human
Gene ID (Human)	3110
Gene ID	3110
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 250-350 of human MNX1/HB9/HLXB9 (P50219).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P50219
Protein Weight	41kDa
Shipping	Ice bag