

Moc3 Rabbit Pab - A0A0N7Kjv6

Item Number: CM0024069

Introduction

MOC3 Rabbit pAb, a polyclonal antibody targeting Oryza sativa MOC3 (SWISS A0A0N7KJV6), validated for WB and ELISA applications in rice research.

[Learn More](#)

Parameter	Value
Product Name	MOC3 Rabbit pAb
Remarks/Alias	WUS; OsWUS; WOX1B; OsWOX1B; Osj_015793; MOC3
Species	Oryza sativa
Gene ID	107275471
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 98-289 of Oryza sativa MOC3. (Q7XM13)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	A0A0N7KJV6
Protein Weight	31kDa
Shipping	Ice bag