


Mocs3 Rabbit Polyclonal Antibody For Wb And Elisa - 095396

Item Number: CM0027988



Introduction

High-quality rabbit polyclonal antibody targeting human MOCS3, validated for WB and ELISA. Recognizes human and mouse, suitable for research on tRNA thiolation and molybdenum cofactor biosynthesis.

[Learn More](#)

Parameter	Value
Product Name	MOCS3 Rabbit pAb
Remarks/Alias	UBA4; MOCS3
Species	Human
GeneID	27304
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 271-460 of human MOCS3 (NP_055299.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95396
Protein Weight	50kDa
Shipping	Ice bag