

Mon1A Rabbit Polyclonal Antibody For Wb And Elisa - Q86Vx9

Item Number: CM0022175



Introduction

MON1A Rabbit Polyclonal Antibody for WB and ELISA. Reacts with Human MON1A, a vacuolar fusion protein involved in membrane trafficking and RAB7 activation. Q86VX9.

[Learn More](#)

Parameter	Value
Product Name	MON1A Rabbit pAb
Remarks/Alias	SAND1; MON1A
Species	Human
GeneID (Human)	84315
GeneID	84315
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human MON1A (Q86VX9)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q86VX9
Protein Weight	73kDa
Shipping	Ice bag