

Monkeypox Virus (Mpxv) E8L Rabbit Mab - Q8V4Y0

Item Number: CM0004712

Introduction

Monkeypox virus (MPXV) E8L rabbit monoclonal antibody validated for WB and ELISA. Recognizes the cell surface-binding protein OPG105 (Q8V4Y0), critical for virion attachment. Ideal for monkeypox virus research.

[Learn More](#)

Parameter	Value
Product Name	Monkeypox virus (MPXV) E8L Rabbit mAb
Species	Monkeypox virus
GeneID	929017
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-275 of human Monkeypox virus (MPXV) E8L. (NP_536532.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Monkeypox virus
SWISS	Q8V4Y0
Protein Weight	35kDa
Shipping	Ice bag