

Monoclonal Antibody Panel For Human And Cynomolgus Monkey Treg Detection By Flow Cytometry - Human/Monkey Treg Panel

Item Number: CM0012219

Introduction

Monoclonal antibody panel designed for flow cytometric identification and analysis of regulatory T cells (Tregs) in human and cynomolgus monkey samples. Ideal for immunophenotyping and immune monitoring research.

[Learn More](#)

Parameter	Specification
Product Name	Human/Monkey Treg Panel
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
Shipping	Ice bag