

Monomethyl-Histone H3 (K14) Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, If, Elisa

Item Number: CM0000369



Introduction

Rabbit polyclonal antibody specific for monomethylated histone H3 at lysine 14. Recognizes human, mouse, and rat H3.1 and H3.1t variants. Validated for Western blot, immunofluorescence, and ELISA. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-K14 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; MonoMethyl-Histone H3-K14
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide; A synthetic monomethylated peptide around K14 of human MonoMethyl-Histone H3-K14 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag